

# Magnetic bead purification of antibodies from cell broth at pilot-scale

Nils A. Brechmann<sup>a,b</sup>, Per-Olov Eriksson<sup>c</sup>, Kristofer Eriksson<sup>a,d</sup>, Sven Oscarsson<sup>e</sup>, Atefeh Shokri<sup>a,b</sup>, Göran Hjälmar<sup>a,d\*</sup>, Véronique Chotteau<sup>a,b\*</sup>

**AdBIOPRO**

Centre for Advanced BioProduction by Continuous Processing

**lab on a bead**  
Beadsurface Biotechnology

## AIM OF THE STUDY

The aim of this study was to develop an efficient mAb capture step in cell broth using magnetic beads, showing a proof of concept of suitability for industrial manufacturing.

## INTRODUCTION


The increasing demand of monoclonal antibodies (mAbs) has placed a significant burden on the downstream process of mAb manufacturing. Protein A column chromatography has limitations with adsorption time and an inability to process cell broth.

Magnetic one-step batch separation eliminates centrifugation and filtration, and fuses the steps of clarification, purification, and concentration. This minimizes the process costs while providing high product purity in a single step.

## RESULTS

Fast adsorption rates were obtained with an average of 88% and 97% mAb captured after 30 min and 1 hour, respectively. The presence of cells, with viabilities of 90% (run B1) or 76% (run B2), had no influence on the process.

From the 24.5 g of mAb present in the 16.25 L cell sample (run B2), a total purification yield of 86% was obtained. The mAb was eluted from the magnetic beads efficiently in a volume comparable to column chromatography. Purity and yield were also comparable to traditional protein A column chromatography (not shown).



Magnetic bead adsorption efficiency




Fig. 1 Adsorption efficiency from run B1 (red/solid) with mAb concentration 1.31 g/L (0.8 L beads) and run B2 (blue/dashed) with 1.51 g/L (1 L beads).

SDS-PAGE uptake analysis




Fig. 2 Remaining mAb during adsorption phase from run B2 at different time points (non-reduced SDS-PAGE).

Magnetic bead elution profile




Fig. 3 Elution profile from run B2 as measured at  $A_{280}$ .

## CONCLUSION

- Magnetic protein A beads were used successfully at pilot-scale to purify mAb directly from non-clarified cellbroth
- The use of magnetic beads eliminates the need for centrifugation/filtration and shortens process time

Purification efficiency




Fig. 4 Purification efficiency showing adsorption and yield percentage.

## Acknowledgements

This work was supported by VINNOVA [grant number 2016-02398] and AdBIOPRO, funded by VINNOVA [grant number 2016-05181]. We also thank Selexis, Switzerland, for the model cell line and Irving Scientific, USA, for the base culture medium.

<sup>a</sup> AdBIOPRO, VINNOVA Competence Centre for Advanced Bioproduction by Continuous Processing, Sweden

<sup>b</sup> Cell Technology Group (CETEG), Dept. of Industrial Biotechnology, School of Engineering Sciences in Chemistry, Biotechnology and Health, KTH-Royal Institute of Technology, SE-10691, Stockholm, Sweden

<sup>c</sup> PE Bioprocess Consulting AB, Strängnäs, Sweden

<sup>d</sup> Lab on a Bead AB, Virdings allé 28, SE-754 50 Uppsala, Sweden

<sup>e</sup> Department of Organic Chemistry, Stockholm University, Stockholm, Sweden

\*Corresponding authors: veronique.chotteau@biotech.kth.se, goran.hjalmar@labonabead.se